



### lonplusឹ

Based on more than 30 years of research and experience, lonplus develops and builds innovative instruments for radiocarbon sample preparation and Accelerator Mass Spectrometry (AMS). Our instruments are highly automated and provide excellent reproducibility and stability that allows our customers to deliver outstanding results in their <sup>14</sup>C and AMS applications. Versatility and a user-friendly design are achieved through excellent engineering.

Ionplus stands for high-quality scientific instruments made in Switzerland as well as excellent customer service. We are dedicated to providing our customers all over the world with the best solutions for <sup>14</sup>C analysis and cutting-edge AMS technology for a wide range of applications. Ionplus offers virtually the entire range of dedicated AMS laboratory and measurement equipment from one source: AMS instruments, fully automated graphitization systems, gas interface systems, automated carbonate handling systems, pneumatic sample presses, vacuum lines for sealing tubes and a range of accessories for all products.

As of 2021, lonplus introduces the latest innovation in AMS: our new low-energy AMS system LEA. It runs at just 50 kV accelerator voltage and is the most compact AMS in the world with dimensions of just 2.6 m x 1.8 m (footprint:  $4.7 \text{ m}^2$ ).



### Scientific Instruments for Radiocarbon Dating and Accelerator Mass Spectrometry

### THE BEST SOLUTIONS FOR ALL AMS APPLICATIONS.

Find the best solutions for your specific needs in  $^{14}\mathrm{C}$  and AMS analysis.

### ARCHAEOLOGY High precision



Archeological samples together with samples from paleoclimatology and radiocarbon calibration projects require highest precision, reproducibility and reliability.

### ENVIRONMENTAL SCIENCE Small sample sizes



From carbon cycle studies to biogeochemistry, lonplus offers the right equipment to deal with challenging sample sizes and sample matrices from a wide variety of research areas.

### MARINE RESEARCH Handling of carbonates and DIC



From sediments to water samples, from the analysis of single foraminifera to dissolved inorganic carbon, lonplus offers equipment for fully automated sample treatment, acidification and sampling.

### MATERIALS SCIENCE Simultaneous $\delta^{13}C$ measurements



As certification of biofuels and plant-derived materials becomes more important, we offer a fast, reproducible all-in-one solution for simultaneous <sup>14</sup>C and δ<sup>13</sup>C analysis. We are there for you – from consulting on application requirements and product combinations to trainings and support in your daily work.

### BIOMEDICINE High throughput, carbon quantification



Ionplus combines the excellent sensitivity of AMS with the ultrafast and automated <sup>14</sup>C analysis required for highthroughput microtracing studies.

### FORENSICS High precision, simultaneous $\delta^{13}$ C measurements



Forensic applications span a wide variety of sample backgrounds and typically rely on the radiocarbon "bomb peak", which can deliver valuable information, in particular when combined with stable isotope information.

### **GEOLOGY** <sup>14</sup>C, <sup>10</sup>Be and <sup>26</sup>AI measurements



Burial and in-situ exposure dating, glacial erosion studies or the determination of denudation rates are just a few of the many applications accessible with the MILEA/MILEA light system through the measurement of <sup>10</sup>Be, <sup>14</sup>C and <sup>26</sup>Al.

### ANTHROPOGENIC NUCLIDES Tracing and monitoring



Tracing and monitoring of anthropogenic I, U and Pu nuclides is a steadily growing field. MILEA offers the highest performance for these applications at lowest energies. For recommended product combinations for most diverse AMS applications please contact us or visit www.ionplus.ch

## 

All lonplus products can be either used as stand-alone instruments or in combination with each other to streamline processes in your <sup>14</sup>C or AMS laboratory. The combination of lonplus instruments ensures the best control

MILEA

Multi-Isotope

MILEA light Multi-Isotope

MICADAS Mini Carbon

Dating System

Low Energy AMS

Accelerator Control

Lab Management

Automated Magazine

LEA

ACS

LMP

Package

AMC

Changer

**INTERFACES** Versatile Gas Measurements

Software

onplus

Low-Energy AMS light

Low-Energy AMS

over variability and allows you to obtain precise and reproducible outcomes. Take full advantage of our integrated Lab Management Package LMP by combining lonplus instruments: virtually all steps of the sample processing can be recorded in a database for convenient and safe logging.

	10 - 12	26	<b>-</b> <b>HTI</b> High Through Interface
•	13	27	GIS Gas Interface System
•	<b>-</b> 14 - 16	28	AGE 3 Automated G Equipment
		29	CHS 2 Carbonate Handling Syst
-•	20 - 22	30	<b>IRMS</b> Isotope Ratio Spectrometer
	23	31	PSP Pneumatic Sample Press
	24	32	<b>FED</b> Ferrum Disper
-•	25	33	<b>TSE</b> Tube Sealing Equipment

nput

Graphitization

tem

Mass

enser

















### AMS Øverview

Comparison of Ionplus AMS Instruments



### - MILEA light





### WIRE DRAWING

NOMINAL VOLTAGE

ISOTOPES

INTENDED USE

PRECISION ON <sup>14</sup>C/<sup>12</sup>C

BLANK ON <sup>14</sup>C/<sup>12</sup>C

DIMENSIONS

NOMINAL WEIGHT

TYPICAL POWER CONSUMPTION

COOLING

### MILEA

300 kV
<sup>10</sup> Be, <sup>14</sup> C, <sup>26</sup> Al, <sup>41</sup> Ca, <sup>129</sup> I U, Pu and
other actinides
Applications of virtually all AMS
isotopes with a focus on $^{\mbox{\tiny 129}}\mbox{I}$ and
actinides
2 ‰
<2E-15
3.5 m x 7 m
(24.5 m <sup>2</sup> )
15 000 kg
Max. 60 kW (Isotope-dependent
Water cooling, 50 kW

### MILEA light

300 kV
$^{10}\text{Be},^{14}\text{C},^{26}\text{Al},^{41}\text{Ca}$ available as an
option
Applications in geology and earth
system science
2 ‰
<2E-15
3.0 m x 4.8 m
(14.4 m <sup>2</sup> )
12 000 kg
Max. 60 kW (Isotope-dependent)
Water cooling, 50 kW

### MICADAS

200 kV
<sup>14</sup> C
All <sup>14</sup> C applications, including
ultra-high precision (1 ‰)
2 ‰ (1 ‰ possible)
<2E-15
3.2 m x 2.6 m
(8.3 m <sup>2</sup> )
4 500 kg
<2.5 kW
Air



### LEA

50 kV
<sup>14</sup> C
All $^{\rm 14}\rm C$ applications above 2 $\%$ precision
and dedicated gas (biomedical, source
apportionment, etc.)
2 ‰
ca. 2E-15
2.6 m x 1.8 m
(4.7 m <sup>2</sup> )
3 000 kg
<2.5 kW
Air

### MILEA Multi-Isotope Low-Energy AMS

The world's most innovative multi-isotope AMS system



### **KEY FEATURES**

- Isotopes: <sup>10</sup>Be, <sup>14</sup>C, <sup>26</sup>Al, <sup>41</sup>Ca, <sup>129</sup>U, Pu and other actinides
- Fast switching between different isotopes
- Optimized ion optics for all isotopes
- Ion source can be equipped with a gas interface for radiocarbon samples
- Fully automated gas measurements with GIS
- Fast magazine changes for continuous measurements without breaking the vacuum or cooling down any part of the ion source
- Low space requirements for a multi-isotope facility through very compact design
- Vacuum insulated accelerator terminal without any moving parts, no SF<sub>6</sub> needed
- Minimal maintenance

In a collaboration, lonplus and ETH Zurich have developed a next generation multi-isotope AMS system at low energies: MILEA. Covering <sup>10</sup>Be, <sup>14</sup>C, <sup>26</sup>Al, <sup>41</sup>Ca, <sup>129</sup>I, U, Pu and other actinides, the new instrument combines the established accelerator and ion source technology of MICADAS with the well-proven concept of the high-energy spectrometer layout of the ETH "TANDY" instrument.

With a footprint of just 3.5 m x 7 m, the space requirements of this new AMS instrument are very low. The accelerator is based on the vacuum insulated MICADAS design, which has been upgraded to support up to 300 kV. The low energy spectrometer of the new instrument comprises an achromatic combination of a 90° electrostatic and magnetic deflector. The layout of the high energy side is inspired by the "TANDY" spectrometer (90° magnetic, 120° electrostatic and 110° magnetic deflector). A quadrupole triplet after the accelerator unit provides similar ion optical conditions for all measured isotopes and facilitates tuning. A new, improved low-noise  $\Delta E$ - $E_{res}$  gas ionization chamber at the back end of the new multi-isotope AMS system provides outstanding separation and identification of interfering particles.

The MILEA prototype instrument built and tested in 2017 has shown excellent performance in all tested applications, similar to or exceeding the performance of its higher energy predecessor at ETH. In combination with the new Accelerator Control Software ACS, MILEA is currently not only the most compact but also the most user-friendly multi-isotope AMS system in the world.

Contact us to learn more about the exciting possibilities with MILEA.











### MILEA Multi-Isotope Low-Energy AMS

### MILEA light Multi-Isotope Low-Energy AMS light





### SPECIFICATIONS

- Isotopes: <sup>10</sup>Be, <sup>14</sup>C, <sup>26</sup>Al, <sup>41</sup>Ca, <sup>129</sup>I, U, Pu and other actinides
- Hybrid cesium negative sputter ion source for solid samples and CO<sub>2</sub> gas
- Random access sample changer with 40 positions
- 300 kV accelerator on a vacuum insulated high voltage platform with a low maintenance solid-state power supply
- Helium stripping
- High mass selectivity, abundance sensitivity of  $<5\cdot10^{-13}$  for U
- 7 Faraday cups and integrators on the HE side covering the entire range of measured currents (1pA – 300µA)
- High resolution  $\Delta E$ - $E_{res}$  gas ionization chamber with absorber cell and two parameter data acquisition
- Height of beam line: 1.15 m for simple maintenance
- Dimensions and weight: 3.5 m  $\times$  7 m x 2 m, 15'000 kg
- No SF<sub>6</sub> insulation gas needed



Scan the QR code or visit www.ionplus.ch and take a 360° virtual product tour.



MILEA light is designed for highest performance measurements of <sup>10</sup>Be, <sup>14</sup>C and <sup>26</sup>Al. It is the simpler, smaller and more affordable alternative to MILEA and the perfect instrument for your demanding applications in geology and all radiocarbon research. The compact light version is largely based on the MILEA principles without an ESA on the lowenergy side and an optimized high-energy setup with a 90° magnetic, 110° electrostatic and another 110° magnetic deflector.

Contact us for more information on our new MILEA light system.

### SPECIFICATIONS

Same as MILEA with the following differences:

- Isotopes: <sup>10</sup>Be, <sup>14</sup>C, <sup>26</sup>Al, (<sup>41</sup>Ca, optional)
- Abundance sensitivity on U: not applicable
- 5 Faraday cups and integrators on the HE side covering the entire range of measured currents (1pA – 300µA)
- Dimensions and weight: 3.0 m x 4.8 m x 2 m, 12'000 kg









- Based on the established MILEA concepts
- Very compact design for a <sup>10</sup>Be, <sup>14</sup>C, <sup>26</sup>Al instrument
- Optimized ion optics for all/isotopes
- Ion source can be equipped with a gas
- interface for radiocarbon samples

### MICADAS Mini Carpon Dating System

The most compact high-performance



### **KEY FEATURES**

- Simple and fast tuning
- High measurement stability over long time
- Fully automated gas measurements with GIS
- Fast magazine changes for continuous measurements without breaking the vacuum or cooling down any part of the ion source
- Low space requirements through very compact design
- Extremely low power consumption of 2.5 kW
- Fully air-cooled system, no cooling water needed
- Hybrid cesium negative sputter ion source for solid and gas cathodes
- Vacuum insulated accelerator terminal without any moving parts, no SF<sub>6</sub> needed
- Minimal maintenance

Highest performance with the world's smallest AMS system: The Mini Carbon Dating System MICADAS is a true precision instrument for your <sup>14</sup>C applications. With its permanent magnet and new design, MICADAS is also the most energy efficient AMS in the world and has the lowest infrastructure requirements.

The first prototype of MICADAS was developed and built by the Laboratory of Ion Beam Physics at ETH Zurich in 2004. Since then, more than 20 instruments have been built and delivered to customers worldwide. MICADAS has demonstrated highest performance and reliability and has become the new standard in <sup>14</sup>C-AMS.

With its dimensions of only 3.2 m x 2.6 m x 2 m, MICADAS is the most compact commercially available <sup>14</sup>C-AMS system in the world. Its helium stripping offers a very high transmission of up to 47 % and outstanding measurement stability, thus significantly reducing the need for retuning. The MICADAS hybrid cesium sputter ion source is equipped with a random-access sample changer that holds up to 40 graphite or gas cathodes. Stable ion beam currents of 50 to 150  $\mu$ A and 10 to 20  $\mu$ A C<sup>-</sup> are readily achieved in routine operation with solid and gas samples, respectively.

The acceleration potential of 200 kV is provided by a solidstate power supply without any moving parts, the terminal is vacuum insulated – no  $SF_6$  or other insulation gases are required. A state of the art gas ionization detector with low noise and virtually no degradation provides the most reliable detection of <sup>14</sup>C ions. With this configuration, blanks older than 50'000 radiocarbon years are readily obtained.

In conjunction with the Gas Interface System GIS, MICADAS performs fully automated gas measurements with an auto sampler, an Elemental Analyzer or CO<sub>2</sub> filled glass or quartz tubes. MICADAS is therefore also the most powerful choice for your small samples and high throughput applications.











### MICADAS Mini Carpon Dating System









### SPECIFICATIONS

- Helium stripping (up to 47 % <sup>14</sup>C-transmission)
- Negative ion currents of 50 to 150 µA on full-sized graphite samples and 10 to 20 µA on gas samples (10 µg carbon or more)
- 200 kV accelerator on a vacuum insulated high voltage platform with a low maintenance solid-state power supply
- Dating of samples back to more than 50'000 radiocarbon years
- Machine blank up to 68'000 radiocarbon years
- Random access sample changer with 40 positions
- Dimensions and weight: 3.2 m x 2.6 m x 2.2 m, 4'500 kg
- 2.5 kW average power consumption
- No cooling water or  $SF_6$  needed



Scan the QR code or visit www.ionplus.ch and take a 360° virtual product tour.



lonplus is dedicated to manufacturing instruments of excellent quality and best performance to enable the exciting science of our customers. Our ambitious goal is to expand the field of AMS by making our technology accessible to all users through user friendly instruments and a high degree of automation and standardization.

One of the key aspects for us to achieve this goal is to offer complete packages and turnkey solutions to our customers. An AMS lab can only work if all aspects from sample preparation to measurement, data handling and data evaluation are taken care of. In addition to the complete package of hardware and software, lonplus is offering extensive training sessions and comprehensive factory and site acceptance testing included in the package.

To lower the hurdle for an AMS acquisition even further, we design our systems in the most compact way with great attention given to the lowest possible site requirements and long maintenance intervals, comparable to other traditional mass spectrometers. To round off the lonplus package, our ACS software allows users to work with all our AMS systems in a state-of-the-art manner.

With our brand new low-energy <sup>14</sup>C system LEA, we introduce the next generation AMS system. LEA is the evolution of the well-established concepts of MICADAS. This new system is designed to become the successor of MICADAS and our new back bone for dedicated <sup>14</sup>C AMS systems. In order to even better address the needs of our commercially oriented users, we push this technology further still and are currently testing an even more compact setup with a footprint of just 1.8 x 2 m.









### lonplus<sup>2</sup>

LEA

Low Energy AMS

### **KEY FEATURES**

- Lowest space requirements through very compact design
- Simple and fast tuning
- Controlled by new ACS software for ease of use
- Fully automated gas measurement with GIS
- Fast magazine changes for continuous
- measurement without breaking the vacuum or cooling down any part of the ion source
- Extremely low power consumption of less than 2.5 kW
- Fully air-cooled system, no cooling water needed
- Hybrid cesium negative sputter ion source for solid and gas cathodes
- Vacuum insulated accelerator terminal without any moving parts, no SF<sub>6</sub> needed
- Minimal maintenance

Based on the concepts of the well-established MICADAS, lonplus has developed a new AMS system around the patented helium stripping technology with an accelerator voltage of just 50 kV – LEA. LEA's lower accelerator voltage led to a significant reduction in footprint without compromising its performance. With its dimensions of just 2.6 m x 1.8 m and a height of 1.5 m, LEA is less than half the size of MICADAS. And by that the new, most compact <sup>14</sup>C-AMS system in the world.

By using much of the same concepts and components found in MICADAS, the new LEA system also offers the same features like: a hybrid Cs sputter ion source for measuring gas and solid samples, random access samples changer with air lock, allowing continuous measurement and the maintenance-free, vacuum insulated tandem accelerator. The helium stripping at low voltages offers an even higher transmission of over 50 % and unparalleled measurement stability.



The world's newest, most compact AMS instrument

Just like its predecessor MICADAS, LEA is compatible with all other lonplus products, such as the Gas Interface System GIS to perform fully automated gas measurements or the Lab Management Package for data collection, reduction and evaluation.

Ionplus' new Accelerator Control Software (ACS) allows you to take full advantage of the ease of use of LEA: From fast autotuning and task engine based work flows to new maintenance and diagnostics features – running an AMS instrument has never been easier.

### SPECIFICATIONS

- Helium stripping, more than 50 % <sup>14</sup>C-transmission
- Negative ion currents of 50 to 150 uA on full-sized graphite samples and 10 to 20 uA on gas samples (10 ug carbon or more)
- 50 kV accelerator on a vacuum insulated high voltage platform with a low-maintenance solid-state power supply
- Precision with a modern reference material:
   0.2 % or better for <sup>14</sup>C/<sup>12</sup>C and <sup>13</sup>C/<sup>12</sup>C ratios
- Dating of samples back to more than 50'000 radiocarbon years
- Random access samples changer with 40 positions
- Dimensions and weight: 2.6 m x 1.8 m x 1.5 m, 3'000 kg
- 2.5 kW average power consumption
- No cooling water or SF<sub>6</sub> needed









ACS Accelerator Control Software



### TUNING CARD

Access all devices relevant to the ion beam. Tune manually or select auto-tuning for the fastest tuning of a device. Save tuning sets and load settings from previous measurements.

### CHAT & LOG, NOTIFICATIONS

Get an overview of current and past events. Comment on tasks and notify other users.

> . .

> > OUT & VENT

""C
 "

: 🗌 Live

### SELECT YOUR VIEW

Each view contains a number of cards with specific information and controls to help you run your tuning, visualize data, do maintenance or change settings.

### TASK LIST

Schedule tasks such as wake up, tuning, measurement and sleep. Let a wizard task guide you through maintenance. A task planner allows you to save your favorite autotuning and measurement routines.

	Š	Tuning						:	Ion Beam	:	Vacuum		Sample
✓ Successful 2021-11-12 13:23	1	Device	lominal				Actual 1	Actual 2	<sup>™</sup> C	165 Hz	Sensor	mbar	V03
2072 - Auto Tuning - ANA Steerer V	~	Source Cathode Potential	1111		2		-5.428 kV	-4.07 mA	0	200	P01: Fore-Vac. T1-T2 P02: Fore-Vac. T3-T8	5.34e-1	0
✓ Successful 2021-11-12 13:23		Source Ion Potential	32456				-29.714 kV	-2.47 mA	°C-LE	62.8 µA	P05: Lock	8.01e-7	~
2073 - Auto Tuning - LE Magnet Successful 2021-11-12 13:27	~	Source Ionizer Heating	1122			1	13.16 A	16.44 V	2C	200 30.52 μΑ	P06: Source P09: LE Accelerator	2.01e-7 2.36e-7	
		Cesium Reservoir	153.0 °C				153.0 °C	±1	0	100	P10: HE Analyser	1.85e-7	
2078 - Stop ✓ Successful 2021-11-12 13:27		Source Extractor Potential	374		ţ		-2.747 kV	-0.92 mA	°C	0.3397 µA	P12: Detector Press.	1:52e+1	
2083 - Single Run	~	Source Steerer X	2376				29.0 V	53		Ť			Change 1 3
✓ Successful 2021-11-12 14:42		Source Steerer Y	3404	_			41.5 V	73	CH	1.332 nA			Channe
2075 - Auto Measurement - DEM058 Successful 2021-11-12 16:05	~	Source Box-Lens	17562				-16.079 kV	-1.340 mA	0	10			Change
2076 - Sleep		LE Magnet	0			1	0.00 A	0.00 V					
V Successful 2021-11-12 16:06		Pulse A (12C)	7731			1	4.718 kV						
2084 - Wake Up - Immediately		Pulse B ( <sup>13</sup> C)	10493		1		6.405 kV	-	2072 - Auto T	uning - LE Magnet	< 3 of 3		
		Pulse Offset	2207	1			1.346 kV	-	< > 20	21-11-12 13:27	2 2000 0		
▲ COLLAPSE HISTORY	ADD TASK	ANA Steerer Y	4475	_			54.6 V	=:					
2085 - Auto Measurement - DEMO58	~	Accelerator	186.00 kV			3	185.93 kV	0.0 μΑ	0.999	********	• •		-1,
Cycle time Cycle		Stripper Control	27767				0.64 MPa	-	0.998-	•		0	-1.
17 of 20 s 3 of 5 Run Target		HE Magnet	-1000		1		-0.31 A	-0.10 V	an le	0		•	
AC1636731346 42.4711.114 Target code Target type		HE ESA	24702			1	37.69 kV	-37.69 kV	tie opr	******	*******		-B
0x20FEF0C3 Inz		"C						0.3397 uA	N N N N N N N N N N N N N N N N N N N	0			
1 of 12 2 of 5		0						-	0.995-	0			-1.0
2 2 s	je l'								0.994- 📀				•
Recalibrate magazine after pass Eject magazine No No									0.993				.0
Created by demo user									-3000	-2000 -100	0 0 10 Nominal value	00 2000	3000





Each user has their own account and only one user at a time is 'in control'.

### STATUS CARDS

See the most important instrument information such as currents, pressures and device states at a glance.

### PLOT CARD

Visualize data and look at tunings. There is also a corresponding view that allows you to look at graphs and tables of all runs.

### Accelerator Control Software





ACS



### **KEY FEATURES**

- Easiest operation of any AMS system
- Highest degree of automation for tuning and measurement
- Intuitive user concept
- Web GUI Client access through web browser from any device
- Secure logging of users and events
- Logging and monitoring of all readbacks over time
- User chat, built-in lab journal

22

Ionplus' Accelerator Control Software ACS is revolutionizing the operation of AMS systems.

Start any of the Ionplus AMS systems with a scheduled wake-up task, let the instrument tune itself and run batches of samples in the most automated way. Workflows are easily saved to provide the user with the easiest, most reproducible AMS experience. Maintenance work is made easy through maintenance tasks, requiring minimum user interaction to shut down, vent or restart the system while guiding the user through the steps.

ACS runs on the control computers of all our AMS systems and is accessed through a web GUI from any client with a common web browser. The sophisticated user management ensures secure and simple instrument control while logging all user actions and events.

### ACS ELEMENTS

- With the new ACS task engine, all lonplus AMS systems are easily operated through predefined workflows.
   Running an AMS instrument has never been easier.
- All system parameters are continuously recorded to monitor the system status. Every change to the system and every user interaction is logged to allow the best data quality and security.
- ACS and the Lab Management Package LMP access the same database so that sample information and measurement data are easily traced throughout the entire sample processing chain and measurement.
- Maintenance tasks automate the shut-down and start-up steps and guide the user through the manual steps of the most common maintenance processes.

lonplus offers a comprehensive package of hardware and software for the efficient handling of samples, sample preparation information, instrument control, measurements and <sup>14</sup>C data reduction. Our Lab Management Package LMP enhances throughput, reliability and quality management of your <sup>14</sup>C-AMS laboratory.

Detailed information on every sample is recorded with the user's sample information and measurement data. All this data is safely stored in a database, allowing fast access without the risks associated with file storage.

While the LMP is most useful in conjunction with our AMS instruments, the sample management part of it can be operated with AGE 3 instruments or as standalone system e.g. for the management of customer data, sample information and preparation details or for data reduction. Contact us for a system tailored to your specific needs.

### LMP ELEMENTS

- Database and software suite for customer and sample information, lab processing information as well as measurement data
- Label printer and bar-code reader for fast and efficient labeling of samples and logging of sample preparation steps
- LMP and the Accelerator Control Software ACS access the same database so that sample information and measurement data are easily traced throughout the entire sample processing chain and measurement.
- LMP works best with the data reduction tool BATS. BATS is a fast and reliable data reduction tool that allows the user to check and visualize the measurement data in real time and performs calibration of radiocarbon dates.







- Ensures data integrity
- Central data storage and backup
- Lab-wide data sharing e.g. between sample preparation on AGE 3 and <sup>14</sup>C measurement on MICADAS
- Convenient user interface for standardized
   data input





### **KEY FEATURES**

- Fully automated magazine handling of virtually unlimited numbers of samples
- Enables an endless and continuous measurement
- Hot swap of magazines
- Individually programmable safety zone respecting the given environment
- Designed and certified for collaborative and safe operation
- Manual operation possible at any time
- Available as an upgrade to existing Ionplus AMS systems
- Can easy be removed for maintenance work
- Stack of 9 magazines per pallet
- Magazines can be changed individually or complete pallets
- Low maintenance

24

The Automated Magazine Changer AMC was developed to increase sample throughput while further reducing user interaction and idle time for all our AMS systems. In combination with the standard lonplus sample changers, AMC enables continuous measurements on all our AMS systems for the first time.

The AMC setup consists of a collaborative robot and a magazine storage platform which allows the fully automated feeding of our AMS systems with up to 10 magazines (40 sample positions each). Single magazines can be placed individually in the storage or the user may swap the entire platform. All measurement and batch related information is stored in the LMP system and available in the ACS software for task planning.



GIS and its peripherals (Elemental Analyzer EA, CHS 2, IRMS) make the lonplus AMS systems the most versatile <sup>14</sup>C instruments on the market. A few examples of the coupling and measurement possibilities:

 Measurement of virtually any volatile organic gas sample in glass ampoules

GIS' built-in ampoule cracker not only allows the direct measurement of  $CO_2$  but also enables the measurement of freons, CO and a range of other compounds from sealed glass ampoules with GIS + our AMS systems.

• Dating of small carbonate samples

Combine CHS 2 + GIS + any lonplus AMS instrument for fully automated sample preparation (leaching, flushing, acidification) and measurement of small carbonate samples.

 Ultrafast screening of solid/liquid samples (biomedical/materials samples)

For the fastest <sup>14</sup>C analysis in just 8 minutes per sample, run EA + GIS + an lonplus AMS instrument in combination. Use the same coupling for more precise measurements by simply increasing the measurement time. To gain precise  $\delta^{13}$ C or  $\delta^{15}$ N information at the same time, add an IRMS instrument for online sample analysis (EA + IRMS + GIS + any lonplus AMS system). Use the new HTI (with EA) for the fastest analysis time of discrete samples in under 5 minutes per sample.

• Coupling to third-party CO<sub>2</sub>-producing instruments

 $\label{eq:GIS} \begin{array}{l} {\sf GIS} + {\sf any \ lonplus \ AMS \ system \ are \ compatible \ with \ any \ CO_2 \mbox{-} producing \ instrument \ with \ just \ a \ few \ adaptations. \end{array}$ 









HT

### High Throughput Interface

### GIS Gas Interface System





### KEY FEATURES

- Direct measurements of CO<sub>2</sub> conjunction with all lonplus AMS systems
- For highly standardized applications (biomedical, biogenic carbon measurements)
- Highest throughput of ca. 300 samples per day
- Highest automation in conjunction with our new AMC system
- No moving parts, maintenance-free operation

The newly developed HTI is similar to the Gas Interface System (GIS) but sets itself apart with a few distinct features: We have minimized cross-talk between samples, increased sample throughput and streamlined the operation through software. As a result, HTI is the best gas interface for applications such as screening, materials testing (biogenig vs. fossil) and biomedical applications with the need for throughput, traceability and reliability. The current implementation allows for a streamlined coupling of an Elemental Analyzer (EA) to the HTI, automated tuning on gas and simple operation through lonplus' ACS software.

Contact us to find the best solution for your application.

### SPECIFICATIONS

- CO<sub>2</sub> sample sizes of 20  $\mu$ g C 2000  $\mu$ g C
- Highly standardized and automated coupling with an Elemental Analyzer EA
- Fastest analysis: less than 5 minutes per sample
- Cross-talk: ca. 0.2% @ 100 μg C
- 4 auxiliary gas inlets for reference gases



The lonplus Gas Interface System GIS is the most versatile gas handling system for <sup>14</sup>C-AMS measurements of CO<sub>2</sub>. Direct measurements of CO<sub>2</sub> are performed on ultra-small samples of 3 to 100 µg carbon with the GIS + lonplus AMS coupling. Gas measurements are the ideal solution not only for small samples but also for all lower precision samples in screening and high throughput studies. Sample CO<sub>2</sub> is mixed with helium and the mixture is continuously fed into the ion source of any lonplus AMS instrument. All functionalities of the instrument are software-controlled and fully automated for gas measurements without user interaction for 8 to 40 samples. Moreover, the coupling of virtually any CO<sub>2</sub>-producing device is possible through the integrated zeolite trap.

### SPECIFICATIONS

- Handles  $CO_2$  sample sizes between 3 and 100  $\mu$ g carbon
- Versatile couplings with an Elemental Analyzer EA, a carbonate system CHS 2 and an automated tube cracker (tube dimensions: 4.0 mm O.D., length: 70–80 mm)
- Online stable isotope information is obtained through a GIS + IRMS coupling
- Fully automated measurements for up to 8 sealed tube samples and up to 40 samples with EA or CHS 2
- 60–150 samples can be handled per day
- 4 auxiliary gas inlets for reference gases









- Direct measurements of CO<sub>2</sub> conjunction with all lonplus AMS systems
- Fully automated sample handling
- Highest versatility through diverse CO<sub>2</sub> sources
- Highest throughput for <sup>14</sup>C measurements
- Blank and reference gases are conveniently measured from pre-mixed gas bottles

### AGE 3 Automated Graphitization Equipment

CHS 2 Carbonate Handling System





### **KEY FEATURES**

- Sample combustion and graphitization combined in one compact system
- Fully automated no user input required after loading samples
- User-friendly software
- No liquid nitrogen required
- Fast graphitization reaction 120 minutes
- High throughput 21 samples per day

The third generation of the Automated Graphitization Equipment AGE 3 is the most compact graphitization system on the market. Used in over 30 laboratories around the world, AGE 3 combines sample combustion and graphite production for AMS in a fast and efficient way. Organics are combusted with an Elemental Analyzer EA, carbonates are hydrolyzed and sampled with the fully automated Carbonate Handling System CHS 2. AGE 3 systems run completely unattended and deliver excellent repeatability due to a high degree of automation. This also shows in the AMS data, where good repeatability translates into higher precision.

### SPECIFICATIONS

- Required carbon content for regular samples: 1-2 mg
- Required carbon content for small samples: > 0.2 mg
- Samples of up to 200 mg containing > 3 % carbon can be processed
- Produced graphite: 0.2–1.0 mg carbon on 3–5 mg iron
- Processing blank: < 0.002 F<sup>14</sup>C
   (> 50'000 radiocarbon years)
- Cross-talk: < 1 ‰ on 1 mg carbon



CHS 2 is the second generation of our head-space sampling system for carbonates, DIC and liquid combustion samples. It is designed for efficient flushing, oxidizing/hydrolyzing and sampling from septum-sealed vials. CHS 2 combines a heater block for up to 64 samples, an adapted auto sampler, a water trap and an adjustable flow regulator. The system is fully implemented in both the AGE 3 and GIS software. As a new feature, CHS 2 comes with a completely automated tool change between acid syringes and sampling needles. An optional water kit for DIC samples of up to 100 ml is also available. With the additional acid containers, an automated leaching, flushing and sampling is now possible without user intervention.

### SPECIFICATIONS

- Dimensions: 830 x 390 x 650 mm, 35 kg
- Tray 1: 64 sample vials 4.5/12 ml
- Tray 2: 9 sample vials 100 ml
- Adjustable flow 0–300 ml/min
- Adjustable temperature: room temperature to 100 °C









### **KEY FEATURES**

- Fully automated flushing/acid addition/sampling (including tool changes) with AGE 3/GIS
   Heater block (20–100 °C) for up/to 64 Labco Exetainers® or 9 100 ml serum bottles
- Heated acid containers for phosphoric acid, leaching acid and oxidant (20/40 ml borosilicate vials)
- Integrated flow regulator with LCD readout and
   flow-alarm
- Integrated water trap with increased capacity for water samples

29

### IRMS Isotope Ratio Mass Spectrometer

### **PSP** Rneumatic Sample Press





### High-precision $\delta^{13}$ C and $\delta^{15}$ N\*\* values are conveniently obtained during graphitization with an AGE 3 instrument or during gas measurements with GIS. A newly implemented Elementar precisION® IRMS instrument is coupled to our AGE 3 or GIS system in order to acquire precise and accurate stable isotope information online. Gain new insight into your samples with this convenient coupling. Applications range from archaeological or forensic samples to materials testing and tracer studies.

- \* The IRMS instrument is a third-party instrument, interfacing with lonplus AGE 3 and GIS instruments.
- \*\*  $\delta^{15}N$  values are obtained in conjunction with an Elemental Analyzer only.

### SPECIFICATIONS

- Dimensions: 595 x 460 x 650 mm, 102 kg
- Typical mass range at 3 kV: 1-76 amu
- Typical split: 3 % for IRMS, remainder for AGE 3/GIS
- Mass resolution:  $> 110 \text{ m}/\Delta \text{m} @ 10 \%$  valley separation



Graphite cathodes are pressed reliably, reproducibly and conveniently with the Pneumatic Sample Press PSP. By the push of a button, an easy to clean pin presses samples into the back of cathodes. By pressing cathodes from the back, surface contamination is significantly reduced and reproducibility of sample currents is improved due to a well-defined graphite position. Cathode holders for cathodes of all AMS manufacturers are available for PSP.

PSP helps you save time in the preparation of AMS cathodes and plays a key role in making high-precision measurements possible.

### SPECIFICATIONS

- Automated pressing by the push of a button
- Adjustable force of 100–800 N for any carbon/catalyst ratio
- Defined pressing time of 1.5 seconds per sample



### KEY FEATURES

- Performs  $\delta^{13}C$  and  $\delta^{13}N$  measurements during graphitization or gas measurements
- Fully automated tuning
- Full functionality for other isotopes such as  $\delta^{18}$ O,  $\delta^{33}$ S and  $\delta^{34}$ S in conjunction with different combustion/pyrolysis setups







- Reproducible, reliable and fast pressing of graphite and other materials for AMS cathodes
  Samples are pressed from the back, resulting in low surface contamination and reproducible currents in the AMS measurement
- Cathode holders for lonplus, NEC and High
   Voltage cathodes are available

### FED Ferrum Dispenser

### **TSE** Tube Sealing Equipment





The iron dispenser FED dispenses a well-defined and reproducible amount of metal catalyst for AGE 3 systems and other graphitization lines. The fast and easy handling of FED saves time and its reproducibility provides the basis for high-precision <sup>14</sup>C measurements. Manually operated and virtually wear-free, FED requires no additional equipment. Repeatability tests indicate that iron masses of typically 4–5 mg\* are obtained readily and reliably with a variability of  $\pm$  2 %.

\*The dispensed mass depends on the mesh size of the iron powder.

### **SPECIFICATIONS**

- Dispenses 4–5 mg of iron powder with a typical variation of ±2 % for a 325 mesh size
- Dosing of iron in ca. 5 seconds per tube
- Accepts any 8 mm O.D. culture tube

KEY FEATURES
 Reproducible dosing of iron or other metallic catalysts for graphitization reactions

- Manually operated
- Maintenance-free



With its compact design, the Tube Sealing Equipment TSE serves as a simple vacuum line to crack ampoules, split samples into several ampoules or to prepare and seal samples for combustion in quartz tubes. TSE is also a convenient and versatile instrument for <sup>14</sup>C preparation laboratories without their own AMS capability. Gas ampoules prepared by TSE can be stored and shipped for later analysis with GIS and vacuum-sealed graphite cathodes can be shipped and stored indefinitely. TSE is manually operated and equipped with two 9 mm Ultra-Torr<sup>®</sup> ports for sealing of combustion tubes, two 4 mm ports for sealing of quartz and glass tubes and one  $\frac{1}{2}$ " Ultra-Torr<sup>®</sup> port with a bellows tube cracker for 9 mm sample tubes. The vacuum line also comprises two calibrated volumes (corresponding to 700 and 2'100  $\mu$ g carbon as CO<sub>2</sub> at 1 bar), two pressure transducers (0-3'000 mbar) with digital readouts and a Peltier cooler for water removal.

### SPECIFICATIONS

- 2 independent pressure transducers
- Calibrated volumes for 700 and 2'100  $\mu g$  carbon
- Peltier cooler for water removal at -20 °C









- Easy to use vacuum line
- 2 calibrated volumes with pressure readouts
- Peltier cooler for water condensation
- Display of pressures, carbon masses and temperatures

Besides manufacturing high-quality scientific instruments, lonplus offers also a wide range of services in the field of radiocarbon dating and AMS.

### TRAINING

Benefit from our experience and expertise in radiocarbon dating and AMS. In our in-house laboratory, we offer hands-on trainings and courses on:

- Operating and maintenance of our instruments
- Sample cleaning and sample preparation for best results
- <sup>14</sup>C analysis and data interpretation
- Best practices in <sup>14</sup>C analysis





### CUSTOM APPLICATIONS AND DEVELOPMENTS

You have a non-standard application or special requests? Or need a custombuilt system specific to your analytical questions? Tell us about your unique challenges and we can provide you with custom solutions tailored to your specific needs.









### **TECHNICAL SUPPORT**

Our experienced staff provides prompt assistance and support for any question or technical problem. We are here for you and your questions by e-mail, phone or even through our remote access assistance.

### MAINTENANCE PACKAGES

The lonplus maintenance packages help you avoid costly downtime, reduce repair work and keep your systems up to date and running. Choose the maintenance package that best suits your needs.

### SPARE PARTS AND UPGRADES

Our instruments are designed and built to last. We offer upgrades of equipment and spare parts for all our instruments.

Currently available upgrade options for older systems:

- Helium stripping for MICADAS
- 40 position magazine sample changer for MICADAS
- IT and software upgrades for AGE, GIS and all our AMS systems

Request a quotation for your desired upgrade at info@ionplus.ch.

### CONSUMABLES

To keep your research going, we deliver consumables for all our products such as:

- Cathodes for solid and gas samples
- Vials and tubes for AGE 3, CHS 2 and TSE
- Reference and blank gas mixtures for GIS













# HEAD QUARTERS













CONTACT

Ionplus AG Lerzenstrasse 12 8953 Dietikon Switzerland T +41 43 322 31 60 info@ionplus.ch www.ionplus.ch



### SALES NETWORK

Ionplus has an international network of sales representatives with local and market-oriented partner companies. For a trading company in your region, please check www.ionplus.ch or contact us at info@ionplus.ch.